



QTM Example - Export Cleanse data


1. Once the event detection is finished, press play



2. The parameters will be split into cycle. You have different tools to exclude cycles from your analysis




3. You can for example use the scalpel tool to drag a line over cycles you want to exclude




4. You can switch between the different layout to check all types of parameters





5. Once your cycle cleaning is done, click on continue




6. In the Cycles Comparison window, you can compare parameters that share the same name accross different analysis




8. Select the analysis you want to export




9. Select the report type




10. Press Run



Once the status is set to "Completed" you can open the folder where your report was exported



11. You can open the report in Excel



You have access to all your parameters in the various tabs of the Excel file

